

San Francisco VA FACS FUSION Configuration Chart											updated 12/16/2019
BLUE (488nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
	FSC	None	None	FSC							
1	A	710/50	690LP	685-735	PE-Cy5.5	PerCP-Cy5.5	BB700	PerCP-efluor710	PerCP-A	7AAD-B	
2	B	530/30	505LP	515-545	FITC	GFP	Alexa 488	CFSE	BB515	YFP	
	SSC	488/10		483-493	SSC						
RED (640nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
3	A	780/60	750LP	750-810	APC-Cy7	APC-H7	Ghost APC 780	APC-FIRE	BD FVD 780	Alexa 750 (dim?)	Alexa 790 (dim?)
4	B	710/50	690LP	685-735	Alexa 700	APC-R700	BD FVD 700				
5	C	670/30		655-685	APC	Alexa 647	BD FVD 660				
VIOLET* (405 nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
6	A	780/60	750LP	750-810	BV786						
7	B	710/50	690LP	685-735	BV711						
8	C	670/30	635LP	655-685	BV650						
9	D	610/20	600LP	600-620	BV605						
10	E	525/50	505LP	500-550	BV510	Alexa 430	eFluor506	Ghost Violet510	Pacific Orange	BD FVD 510	BV480
11	F	450/50		425-475	BV421	Sytox Blue	Pacific Blue	Alexa 405	Dead Violet	BD FVD 450	BD V450
Yellow Green (561 nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
12	A	780/60	750LP	750-810	PE-Cy7						
13	B	670/30	635LP	655-685	PE-Cy5	7AAD-A	PerCP-B				
14	C	610/20	600LP	600-620	PE-efluor610	Dead Red	PE-Dazzle 594	PE-CF594	BD FVD 620	PE-Texas Red	Propidium Iodide
15	D	586/15		578-593	PE	BD FVD 570	TdTomato				
Try these spectral viewers to help you design your experiment, and determine fluor compatibility with the machine configuration.											
https://www.bdbiosciences.com/en-us/applications/research-applications/multicolor-flow-cytometry/product-selection-tools/spectrum-viewer											
https://www.biolegend.com/en-us/spectra-analyzer											
https://www.thermofisher.com/us/en/home/life-science/cell-analysis/labeling-chemistry/fluorescence-spectraviewer.html											
General guidance rules: If you have more than 1 marker on the same cell, separate them onto different lasers (e.g. CD3 PerCP and Ly6C APC). Do not put an abundant marker onto your brightest											

San Francisco VA FACS ARIAIU Configuration Chart updated CH 10/29/2018

BLUE (488nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
FSC	None	None		FSC							
1	A	780/60	735LP	750-810	PE-Cy7						
2	B	695/40	655LP	675-715	PerCP	PE-Cy5.5***	PE-Cy5***	PerCP-eFluor710	PerCP-Cy5***	7-AAD	PerCP-Cy5.5***
3	C	610/20	595LP	600-620	PE-eFluor610	Pe-Texas Red	Dead Red	PE CF 594	PE-Dazzle 594		
4	D	576/26	556LP	563-589	PE	Propidium Iodide	Alexa Fluor 555				
5	E	530/30	502LP	515-545	FITC	GFP	Alexa Fluor 488	CFSE	ALDH	DEAB	
	F	488/10		483-493	SSC						

RED (633nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
6	A	780/60	735LP	750-810	Alexa Fluor 750	APC-Cy7	APC-H7	Alexa Fluor 790	Ghost APC 780	APC-FIRE	
7	B	660/20		650-670	APC***	Alexa Fluor 647	Alexa Fluor 610				
	C										

VIOLET* (405 nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
8	A	670/30	635LP	655-685	BV650***						
9	B	530/30	505LP	515-545	Pacific Orange**	Alexa Fluor 430	eFluor506	BV510	GHOST VIOLET 510		
10	C	450/50		425-475	Sytox Blue	Pacific Blue**	Alexa Fluor 405	Dead Violet	BV421		
	D										
	E										
	F										
	G										
	H										

UV* (375nm excitation)											
Detector	Band pass Filter	Longpass filter	Effective range	FLUOROCHROMES							
	A	670/30	635LP	655-685							
	B	530/30	505LP	515-545							
11	C	450/50		425-475	DAPI	eFluor 455UV	Hoescht	Alexa Fluor 350			
	D										
	E										
	F										
	G										
	H										

*The violet and UV lasers cannot be used simultaneously.
 ** In the Violet laser section, Pacific Orange and Pacific Blue have the same excitation and emission spectra (ex405/em455, 500, 551)
 *** In the Blue and Red laser sections, Cy5 and Cy5.5 tandem conjugates are poor choices due to spills, crossactivation and instability. Try to avoid PE-Cy5(.5) conjugates
General guidance rules: If you have more than 1 marker on the same cell, separate them onto different lasers (e.g. CD3 PerCP and Ly6C APC). Do not put an abundant